

## PCR – stretch films Machine stretch films, manual stretch films & cover sheets

Ecologically.  
sustainable. good.

### PCR-Machine stretch films

#### LOOP SAVE 30

- >> 300 % elongation
- >> thickness: 12-35 microns

#### LOOP PROFI P 30

- >> 250 % elongation
- >> thickness: 15-35 microns

#### LOOP P II 30

- >> 200 % elongation
- >> thickness: 15-35 microns

#### LOOP P I 30

- >> 150 % elongation
- >> thickness: 15-35 microns

### PCR-manual stretch films

#### LOOP P I 30

- >> thickness: 15-35 microns

#### LOOP SAVE 30/60

- >> thickness: 10-15 microns

### „Blauer Engel“- cover characteristics:

- size: 1.200 x 1.600 mm
- strength: 30 microns
- 250 pieces / Roll,
- perforated
- Naturally cloudy or
- translucent
- 35 rolls / pallet



### Stretch films with PCR- Recyclate content of 30%

We test on site to see what quality you get brings the most advantages. Make an appointment with us!

### POST-CONSUMER-RECYCLING (PCR)

What does PCR consist of?

- >> Plastic waste generated by end consumers (classic household waste) as well as packaging waste from retail, supermarkets and sales.

**Do you have any further questions?  
We would like to help you!**

**We would also be happy make you  
an offer.**

OptiFol GmbH  
Carl-Zeiß-Str. 9  
97424 Schweinfurt  
Telefon: 0 9721 97850-70  
E-Mail: [info@optifol.de](mailto:info@optifol.de)  
Web: [www.optifol.de](http://www.optifol.de)